



701 Hemlock Road
Morgantown PA 19543
Phone: 610.524.5810 Fax: 610.901.3016
info@AzerSci.com www.AzerSci.com

Certificate of Analysis

1/7/2022

Product No.: ES613-32

Description: Reagent Grade Water, 32oz, 12/case

Lot No.: 27430

Date of Manufacturing: 1/6/2022

Expiration Date: 7/6/2023

RMC-DIH2O Lot #: 27430

| PROPERTY | TEST SPECIFICATION | RESULTS |
|---------------------|----------------------------|------------|
| Appearance | Clear, colorless, odorless | Conforms |
| Resistivity | 10 MΩ /cm | 18.03MΩ/cm |
| Silicates (soluble) | <50 µg/L | 3.2 µg/L |
| Sodium ion | <10.0 µg/L | 0.13 µg/L |
| Chloride ion | <10.0 µg/L | 0.16 µg/L |
| TOC | <50 µg/L | 14 µg/L |
| Bacterial Count | < 10 CFU/ml | 0 CFU/ml |

Country of Origin: United States of America

Latex Free: This product contains no latex

Created By:

Fernando Fana
Quality Associate