



# Certificate of Analysis

2/10/2022

701 Hemlock Road  
Morgantown PA 19543  
Phone: 610.524.5810 Fax: 610.901.3016  
info@AzerSci.com www.AzerSci.com

**Product No.:** B-PFNBF-20

**Description:** Prefilled 10% Neutral Buffered Formalin, 20mL, half-filled, Formalock closure, bulk packaging 32/pack, 24 packs/case

**Lot No.:** 27714

**Date of Manufacturing:** 2/7/2022

**Expiration Date:** 1/1/2024

## 10% Formalin Lot #: 2A010

PROPERTY	TEST SPECIFICATION	RESULTS
PH	6.90-7.10 @ 25°C	6.94@ 25°C
Formaldehyde Concentration	9.5-10.5% (v/v)	10.31%(v/v)
Appearance	Clear, colorless, no observable particulate	Conforms

## 10% Formalin Lot #: 2B001

PROPERTY	TEST SPECIFICATION	RESULTS
PH	6.90-7.10 @ 25°C	7.01@ 25°C
Formaldehyde Concentration	9.5-10.5% (v/v)	10.31%(v/v)
Appearance	Clear, colorless, no observable particulate	Conforms

**Country of Origin:** United States of America

**Latex Free:** This product contains no latex

**Created By:**

James E Osborn  
Quality Associate