



701 Hemlock Road
Morgantown PA 19543
Phone: 610.524.5810 Fax: 610.901.3016
info@AzerSci.com www.AzerSci.com

Certificate of Analysis

2/11/2022

Product No.: B-PFNBF-60

Description: Prefilled 10% Neutral Buffered Formalin, Formalock closure, 60mL, half-filled, bulk, 27/pack 324/case

Lot No.: 27790

Date of Manufacturing: 2/7/2022

Expiration Date: 1/1/2024

10% Formalin Lot #: 2A010

| PROPERTY | TEST SPECIFICATION | RESULTS |
|----------------------------|---|-------------|
| PH | 6.90-7.10 @ 25°C | 6.94@ 25°C |
| Formaldehyde Concentration | 9.5-10.5% (v/v) | 10.31%(v/v) |
| Appearance | Clear, colorless, no observable particulate | Conforms |

Country of Origin: United States of America

Latex Free: This product contains no latex

Created By:

James E Osborn
Quality Associate