



# Certificate of Analysis

2/11/2022

701 Hemlock Road  
Morgantown PA 19543  
Phone: 610.524.5810 Fax: 610.901.3016  
info@AzerSci.com www.AzerSci.com

**Product No.:** PFNBF-1.25G

**Description:** Prefilled 10% Neutral Buffered Formalin, 1.25 gallon, half-filled, Each (MOQ-12)

**Lot No.:** 27648

**Date of Manufacturing:** 2/10/2022

**Expiration Date:** 1/1/2024

## 10% Formalin Lot #: 2A011

PROPERTY	TEST SPECIFICATION	RESULTS
PH	6.90-7.10 @ 25°C	7.01 @ 25°C
Formaldehyde Concentration	9.5-10.5% (v/v)	10.16%(v/v)
Appearance	Clear, colorless, no observable particulate	Conforms

## 10% Formalin Lot #: 2B002

PROPERTY	TEST SPECIFICATION	RESULTS
PH	6.90-7.10 @ 25°C	7.02 @ 25°C
Formaldehyde Concentration	9.5-10.5% (v/v)	10.31%(v/v)
Appearance	Clear, colorless, no observable particulate	Conforms

**Country of Origin:** United States of America

**Latex Free:** This product contains no latex

**Created By:**

James E Osborn  
Quality Associate