



Certificate of Analysis

3/9/2022

701 Hemlock Road
Morgantown PA 19543
Phone: 610.524.5810 Fax: 610.901.3016
info@AzerSci.com www.AzerSci.com

Product No.: B-PFNBF-20

Description: Prefilled 10% Neutral Buffered Formalin, 20mL, half-filled, Formalock closure, bulk packaging 32/pack, 24 packs/case

Lot No.: 28055

Date of Manufacturing: 3/3/2022

Expiration Date: 2/1/2024

10% Formalin Lot #: 2B012

PROPERTY	TEST SPECIFICATION	RESULTS
PH	6.90-7.10 @ 25°C	6.98@ 25°C
Formaldehyde Concentration	9.5-10.5% (v/v)	10.31%(v/v)
Appearance	Clear, colorless, no observable particulate	Conforms

10% Formalin Lot #: 2C001

PROPERTY	TEST SPECIFICATION	RESULTS
PH	6.90-7.10 @ 25°C	6.92@ 25°C
Formaldehyde Concentration	9.5-10.5% (v/v)	10.00%(v/v)
Appearance	Clear, colorless, no observable particulate	Conforms

Country of Origin: United States of America

Latex Free: This product contains no latex

Created By:

Fernando Fana
Quality Associate