



701 Hemlock Road  
Morgantown PA 19543  
Phone: 610.524.5810 Fax: 610.901.3016  
info@AzerSci.com www.AzerSci.com

# Certificate of Analysis

3/2/2022

**Product No.:** ES616-4G

**Description:** Deionized Water, 1 gallon, 4/case

**Lot No.:** 27988

**Date of Manufacturing:** 2/23/2022

**Expiration Date:** 8/23/2023

## RMC-DIH2O Lot #: 27988

| PROPERTY            | TEST SPECIFICATION         | RESULTS    |
|---------------------|----------------------------|------------|
| Appearance          | Clear, colorless, odorless | Conforms   |
| Resistivity         | 10 MΩ /cm                  | 18.21MΩ/cm |
| Silicates (soluble) | <50 µg/L                   | 3.2 µg/L   |
| Sodium ion          | <10.0 µg/L                 | 0.13 µg/L  |
| Chloride ion        | <10.0 µg/L                 | 0.16 µg/L  |
| TOC                 | <50 µg/L                   | 14 µg/L    |
| Bacterial Count     | < 10 CFU/ml                | 0 CFU/ml   |

**Country of Origin:** United States of America

**Latex Free:** This product contains no latex

**Created By:**

James E Osborn  
Quality Associate