



701 Hemlock Road  
Morgantown PA 19543  
Phone: 610.524.5810 Fax: 610.901.3016  
info@AzerSci.com www.AzerSci.com

# Certificate of Analysis

4/25/2022

**Product No.:** ES124310L

**Description:** Blood Bank Saline, buffered, 7.0-7.2 pH range, 10L cube with spigot, Each

**Lot No.:** 28820

**Date of Manufacturing:** 4/21/2022

**Expiration Date:** 10/19/2023

## Buffered Saline Lot #: 2D28818

| PROPERTY   | TEST SPECIFICATION               | RESULTS    |
|------------|----------------------------------|------------|
| PH         | 7.0-7.2 @ 25°C                   | 7.11       |
| Assay      | 0.85% to 0.90%                   | 0.89%      |
| Appearance | Clear, Colorless, No Particulate | Conforms   |
| TAMC       | < 10 CFU/ml                      | < 1 CFU/ml |
| TYMC       | < 10 CFU/ml                      | < 1 CFU/ml |

**Country of Origin:** United States of America

**Latex Free:** This product contains no latex

**Created By:**

James E Osborn  
QA Associate