



701 Hemlock Road
Morgantown PA 19543
Phone: 610.524.5810 Fax: 610.901.3016
info@AzerSci.com www.AzerSci.com

Certificate of Analysis

4/13/2022

Product No.: 10-820

Description: Kit, 12 place Prostate Biopsy, 12x20mL vials half filled with 10% formalin per kit, 25 kits per case

Lot No.: 28691

Date of Manufacturing: 4/12/2022

Expiration Date: 3/1/2024

10% Formalin Lot #: 2C017

PROPERTY	TEST SPECIFICATION	RESULTS
PH	6.90-7.10 @ 25°C	7.02@ 25°C
Formaldehyde Concentration	9.5-10.5% (v/v)	10.46%(v/v)
Appearance	Clear, colorless, no observable particulate	Conforms

10% Formalin Lot #: 2D001

PROPERTY	TEST SPECIFICATION	RESULTS
PH	6.90-7.10 @ 25°C	6.97@ 25°C
Formaldehyde Concentration	9.5-10.5% (v/v)	10.46%(v/v)
Appearance	Clear, colorless, no observable particulate	Conforms

Country of Origin: United States of America

Latex Free: This product contains no latex

Created By:

James E Osborn
QA Associate