

## Certificate of Analysis

**Product Name:** Blood Bank Saline – Buffered Solution

**Description:** Isotonic Saline Solution, 0.85% w/v, Buffered

**Product No.:** ES124320L

**Lot. No.:** 21184

**Date of Manufacture:** 07/08/2020

**Expiration Date:** 01/06/2022

**Product Specifications:**

**Ingredients:** Sodium Chloride, USP Grade  
Disodium Phosphate, Food Grade  
Monosodium Phosphate, Food Grade  
Deionized Water, NCCLS Type I Reagent Grade

**Packaging:** Product undergoes a two-phase sub-micron filtration process as well as UV Sterilization before filling process begins.  
Packaged in 20 Liter Cubetainers (Manufactured by Hedwin Corporation.)

**Materials:**

	<u>Part No.</u>	<u>Lot No.</u>	<u>Exp. Date</u>
<b>Container:</b>	CDC6378	05/20/2020	N/A
<b>Cap:</b>	CAP-1 Gal- RMS-38	202004572023	N/A
<b>Spigot:</b>	COM4401-IW	1078477	N/A
<b>Chemical:</b>	ES1243-1000	0G21142	2022-01-06

**QC Testing:**

LOT NUMBER: 0G21142			
ATTRIBUTE	SPECIFICATION	RESULTS	PASS / FAIL
Appearance	Clear, colorless, no observable particulate	Pass	PASS
pH	7.0-7.2	7.08	PASS
Assay	0.85% to 0.90%	0.88%	PASS
Isotonicity (Macro observation)	No hemolysis	Pass	PASS
Isotonicity (Micro observation)	No RBC aggregates	Pass	PASS
TAMC	< 10 CFU/ml	0 CFU/ml	PASS
TYMC	< 10 CFU/ml	0 CFU/ml	PASS

**Country of Origin:** United States of America

**Latex Free:** This product contains no latex

Created By:  
Fernando Fana  
Quality Associate