

Certificate of Analysis

Product Name: Blood Bank Saline – Buffered Solution

Description: Isotonic Saline Solution, 0.85% w/v, Buffered

Product No.: ES124310L

Lot. No.: 22713

Date of Manufacture: 11/10/2020

Expiration Date: 05/09/2022

Product Specifications:

Ingredients: Sodium Chloride, USP Grade
 Disodium Phosphate, Food Grade
 Monosodium Phosphate, Food Grade
 Deionized Water, NCCLS Type I Reagent Grade

Packaging: Product undergoes a two-phase sub-micron filtration process as well as UV Sterilization before filling process begins.
 Packaged in 10 Liter Cubetainers (Manufactured by Hedwin Corporation.)

Materials:

| | <u>Part No.</u> | <u>Lot No.</u> | <u>Exp. Date</u> |
|-------------------|------------------|----------------|------------------|
| Container: | CDC5894 | 10/06/2020 | N/A |
| Cap: | CAP-1 Gal-RMS-38 | 202009586132 | N/A |
| Spigot: | COM4401-IW | 1078478 | N/A |
| Chemical: | ES1243-1000 | 0K22709 | 05/09/2022 |

QC Testing:

| LOT NUMBER: 0K22709 | | | |
|--|---|----------------|--------------------|
| ATTRIBUTE | SPECIFICATION | RESULTS | PASS / FAIL |
| Appearance | Clear, colorless, no observable particulate | Pass | PASS |
| pH | 7.0-7.2 | 7.14 | PASS |
| Assay | 0.85% to 0.90% | 0.90 | PASS |
| Isotonicity (Macro observation) | No hemolysis | Pass | PASS |
| Isotonicity (Micro observation) | No RBC aggregates | Pass | PASS |
| TAMC | < 10 CFU/ml | 0 CFU/ml | PASS |
| TYMC | < 10 CFU/ml | 0 CFU/ml | PASS |

Country of Origin: United States of America

Latex Free: This product contains no latex

Created By:
 James Osborn
 Quality Associate