

## Certificate of Analysis

**Product Name:** Blood Bank Saline – Standard Solution

**Description:** Isotonic Saline Solution, 0.85% w/v, Non-buffered

**Product No.:** ES124420L

**Lot. No.:** 21354

**Date of Manufacture:** 07/21/2020

**Expiration Date:** 01/20/2022

**Product Specifications:**

**Ingredients:** Sodium Chloride, USP Grade  
Deionized Water, NCCLS Type I Reagent Grade

**Packaging:** Product undergoes a two-phase sub-micron filtration process as well as UV Sterilization before filling process begins.

Packaged in 20 Liter Cubetainers (Manufactured by Hedwin Corporation.)

**Materials:**

	<u>Part No.</u>	<u>Lot No.</u>	<u>Exp. Date</u>
<b>Container:</b>	CDC6378	06/23/2020	N/A
<b>Cap:</b>	CAP-1 Gal- RMS-38	202004572028	N/A
<b>Spigot:</b>	COM4401-IW	1078477	N/A
<b>Chemical:</b>	ES1244-1000	0G21353	2022-01-20

**QC Testing:**

LOT NUMBER: 0G21353			
ATTRIBUTE	SPECIFICATION	RESULTS	PASS / FAIL
Appearance	Clear, colorless, no observable particulate	Pass	PASS
pH	6.5 - 7.0	6.6	PASS
Assay	0.85% to 0.90%	0.87%	PASS
Isotonicity (Macro observation)	No hemolysis	Pass	PASS
Isotonicity (Micro observation)	No RBC aggregates	Pass	PASS
TAMC	< 10 CFU/ml	0 CFU/ml	PASS
TYMC	< 10 CFU/ml	0 CFU/ml	PASS

**Country of Origin:** United States of America

**Latex Free:** This product contains no latex

Created By:  
Fernando Fana  
Quality Associate