

Certificate of Analysis

Product Name: Reagent Grade Water, 32oz 12/cs

Description: Reagent grade, CAP Type I, Type II, CLSI, USP, ASTM III Type C

Product No.: ES613-32

Lot. No.

23229

Date of Manufacture:

2021-01-08

Expiration Date:

2022-07-08

QC Testing:

| LOT NUMBER: 23229 | | | |
|----------------------------|----------------------------|----------------|--------------------|
| ATTRIBUTE | SPECIFICATION | RESULTS | PASS / FAIL |
| Appearance | Clear, colorless, odorless | Pass | PASS |
| Resistivity | ≥10 MΩ/cm | 18.37 MΩ/cm | PASS |
| Silicates (soluble) | <50 µg/L | <0.5 µg/L | PASS |
| Sodium ion | < 10.0 µg/L | 0.03 µg/L | PASS |
| Chloride ion | < 10.0 µg/L | 0.05 µg/L | PASS |
| TOC | < 50 µg/L | <6 µg/L | PASS |
| Bacterial count | < 10 CFU/ml | 0 CFU/ml | PASS |

Current process includes a combination of cold process technologies including reverse osmosis, ion exchange resin processes, UV sterilization, and micro-filtration. Analysis data are reported from the most recent scheduled monitoring.

Country of Origin: United States of America

Latex Free: This product contains no latex

Created By:
James Osborn
Quality Associate