

15mL & 50mL Centrifuge Tubes



Conical 15mL and 50mL centrifuge tubes are made of premium polypropylene for clarity and centrifugal strength. Tubes and caps may be frozen at temperatures as low as -86°C. The maximum centrifugal force is 12,000 x g for the 15mL tubes and 15,000 x g for the 50mL tubes. Caps have a flat top surface for easy marking and labeling.

Features:

- Conical Centrifuge Tubes
- High Clarity Polymers
- Printed Graduations & Marking Spot
- RNase & DNase Free, Sterile

Cat. No.	Description	UOM
ES-9151N	15mL Centrifuge Tube, Flat Cap, Racked	500/case
ES-9152N	15mL Centrifuge Tube, Flat Cap, Bulk	500/case
ES-9501N	50mL Centrifuge Tube, Flat Cap, Racked	500/case
ES-9502N	50mL Centrifuge Tube, Flat Cap, Bulk	500/case