

# ALI-Q™2 Repeat Pipettor

Faster and Easier Aliquoting. No Need to Watch the Lines.  
More Accurate and Reproducible Than by Eye.



Set Aliquot Speed  
(VS) model only



Set Aliquot Volume



Press Aliquot Button  
To Dispense Each Aliquot

## REPEAT PIPETTING PLUS CONVENTIONAL PIPETTING WITH A SINGLE 2-IN-1 CONTROLLER.

- Highly reproducible accuracy with the press of one button
- Use any size and any brand of serological pipets
- Plus aspirate and dispense as you're used to
- Ideal for cell culture, dispensing media, filling plates, tubes, cryovials, etc.

You asked and we listened! The ali-Q™2 offers more versatility with features to make your work even easier.

## NEW FEATURES OF THE ALI-Q™ 2 SERIES:

- Speed dial for adjustable aliquoting speed (VS model only)
- On-board calibration software
- Removable kickstand
- Lighter weight

Cat #	Model	Volume Range	Speed Control (Dial)	Speed/Force	On-board Calibration	Removable Kickstand
	ali-Q™2	0.5 - 5.0 mL	No	High (fixed)	Yes	Yes
	ali-Q™2 LS, Low aliquoting speed	0.3 - 3.0 mL	No	Low (fixed)	Yes	Yes
	ali-Q™2 VS, Variable aliquoting speed	0.5 - 5.0 mL	Yes	Custom adjustable (7 speeds)	Yes	Yes

Wobble Nots and AliQ are registered trademarks of Vistalab Technologies Inc.