

# Saliva Specimen Collection Device



**STERILE**

**1 saliva collector funnel, with 10mL transport tube and 1 white cap per set**

The Saliva Collection Kit is designed for saliva sample collection for further transport and testing of saliva specimens. Each kit comes sterile and sealed with one saliva specimen funnel collection device, one 10 mL transport tube, and one white cap per set. Collection tubes are provided empty and ready to be filled with a buffer or collection or preservative solution.



Cat #	Description	UOM
ES45081-10	Saliva Specimen Collection Device, 10mL Transport tube	500 sets/case