

# Disposable Face Shield



Cat #	Description	UOM
ES86433-025	Disposable face shield, 14"W x 9"D antifog shield, unassembled, anti-microbial frame	50/kits per pack

Protect your eyes, mouth, and nose with Azer Scientific's clear disposable face shield. Easy to put on and comfortable to wear, our shield features two adjustable angle settings, antifog cover and anti-microbial frame. While fully protecting the face, the disposable face shield leaves enough room to be comfortably worn even while wearing glasses and a mask.

- 14" wide and 9" deep
- Two adjustable angle settings
  - straight for standard use and bottom angled out for use with a larger mask/respirator
- Antifog Shield
- Anti-microbial frame
- Fully covers face
- Unassembled
- Made in USA