

# Everlast 2D Rocker w/ Flat Mat



**Features:**

- Everlast™ 5 year warranty
- Stainless steel surface
- 16 lbs. load capacity
- Adjustable speed & tilt angle
- Optional stacking platform

With an unprecedented 5 year limited warranty, our Everlast 2D Rocker is built to last. Twenty four hours a day, seven days a week... No problem! The powerful motor and heavy-duty drive system remain reliable and consistent, even under constant usage in a cold room (min. 0°C) or incubator (max. 65°C).

The rocking speed and tilt angle can easily be adjusted to vary between gentle and aggressive agitation. Both the standard platform and the (optional) stacking platform are constructed of stainless steel, preventing against any rusting or corrosion. The Everlast is designed to provide years or uninterrupted service.



Everlast 2D Rocker shown w/ ES59451-STACK

**Technical Data:**

Speed range:	3 to 80 rpm
Tilt angle:	0 to 20° (user adjustable)
Platform size:	16 x 12 in. (40.5 x 30 cm)
Dimensions:	16 x 13 x 6 in. (40.5 x 33 x 15.25 cm.)
Load capacity:	12 lbs. (4.45 kg)
Operating temp:	0°C to +65°C
Weight:	20 lbs. (9 kg)

**Cat. No.**

**Description**

ES59451	Everlast 2D Rocker with flat mat, 115V
ES59451-STACK	Stacking platform with flat mat, 3.5" clearance