

- Complete prostate biopsy collection kit for distribution directly to clinical site
- Packaged ready-to-use with eight, twelve, or sixteen 20 mL Formalock™ containers
- Each container is half-filled with 10% neutral buffered formalin
- Each kit contains an insert with a pre-labeled map of the prostate for assigning biopsy locations
- Collection and transport instructions conveniently printed inside kit lid



Don't see what you are looking for here?
We are able to customize your kit with
components such as biohazard bags and
formalin absorbent pads.

Contact us to learn more!

16 Place Insert



12 Place Insert