

# Low Speed Orbital Shaker



#### Features:

- Low speed, 3 to 70rpm
- Gentle, orbital motion for sensitive gels/blots
- Incubator and cold room safe

The Low Speed Orbital Shaker features a gentle orbital motion with speed control down to as little as 3 rpm, making it the ideal orbital choice for low speed applications, such as blotting, washing, staining/destaining, etc. A large 14 x 12 in. platform is included along with a non-slip rubber mat for holding trays, dishes, plates, flasks and other flat bottomed vessels.

#### Technical Data:

Speed range:	3 to 70rpm
Orbit:	19mm (3/4 in.)
Platform Dimensions:	14 x 12 in. (35 x 30 cm)
Dimensions (WxDxH):	14 x 12 x 8 in. (35 x 30 x 20 cm)
Weight:	12 lbs. (5.5 kg)
Operating temp.:	4°C to 50°C
Electrical:	120V or 230V, 50-60Hz

#### Cat. No.

#### Description

ES59430

Low Speed Orbital Shaker