

Mini Microplate Centrifuge

- Spins down samples in PCR plates
- Accepts all popular PCR plates including non-skirted, semi-skirted, and skirted
- Unique rotor allows user to insert plates into carriers set at a 75 angle before centrifugation. This eliminates the concern of sample loss while inserting plates into centrifuge.
- Once the lid is closed, carriers will reposition into a horizontal position by centrifugal force

Speed: 2,550 rpm (600 xg)

Capacity: 2 x 2 PCR plates
2 x 2 Microtiter plates
24 x 0.2 mL PCR strips

Deceleration: 4 seconds

Dimensions: 9.2 x 10.2 x 7.75 (in)
23 x 26 x 19.7 (cm)

Weight: 9 lbs (4.1 kg)

Electrical: 115 or 230 V, 50-60 Hz, 250W



Before & After
Centrifugation
(angled position)



During
Centrifugation
(horizontal position)

Cat. No.	Description	UOM
ES59770	Mini Microplate Centrifuge	Each
ES59772	Adapter 96 x 0.2 mL/ 12 x PCR strip	Each
ES59774	Rotor for Deep Well Plates	Each