

Mini Centrifuge



Storage Compartment



Includes Both Rotors

Features:

- Instant rotor attachment
- Eight position microbe rotor
- Storage compartment for PCR rotor and 0.5mL & 0.2mL adapters
- Decelerates in just 1 second
- Four position PCR strip rotor

The Mini Centrifuge is a personal centrifuge that virtually fits in the palm of your hand. Ideal for quick spin downs of microtubes and PCR tubes, the Mini Centrifuge is extremely easy to use. Simply close the lid and the rotor instantly reaches 6,000 rpm. Open the lid, and the improved braking system quickly and smoothly decelerates the rotor to a complete stop in just one second.

An eight-position microbe is included and comes installed. Also included is a second rotor with increased capacity for four 0.2mL PCR strips (or 32 x 0.2mL tubes) which is stored in a convenient compartment on the bottom of the centrifuge. This compartment is also used for securing the included tube adapters for 0.5mL tubes.

A new rotor attachment system eliminates the requirement for set screws and tools during rotor exchange. Simply remove and attach rotors by hand, pulling or pressing the rotor onto the motor shaft.

Technical Data:

Speed:	6,000rpm (2,000xg)
Capacity:	8 x 1.5 / 2.0mL PCR tubes 32 x 0.2mL PCR tubes 4 x PCR strips (8 x 0.2mL)
Dimensions:	5.7 x 6.9 x 4.5 in. (14.5 x 17.8 x 11.4 cm)
Weight:	2.2 lbs (1 kg)
Electrical:	100 - 240V, 50-60Hz

Cat. No.

Description

ES89421B	Mini Centrifuge, Blue
ES89421-C	Mini Centrifuge, Clear
ES89421G	Mini Centrifuge, Green
ES89421-PL	Mini Centrifuge, Purple
ES89421R	Mini Centrifuge, Red