

# Mortexer Vortex Mixer w/ Multi-Head™



## Features:

- The only vortexer available with Multi-Head™
- Powerful Motor for INSTANT Vortexing
- Q-Drive™ Dynamic Balance System
- Touch or continuous operation



The Mortexer™ is MORE than your standard vortex mixer. The unique design of the Multi-Head™ (included) features a traditional cup head for general purpose mixing, plus holds up to eight microcentrifuge tubes.

The drive system is balanced with Q-Drive™ technology, providing a smooth vortexing action, while minimizing noise and excessive vibration. This results in quieter, more efficient and longer lasting operation. Increased speed, up to 3400 rpm produces an instant vortex of even the largest sample sizes, including full 50mL tubes.

Other features include: "no-walk" suction feet, cold room/incubator compatibility and spill-proof electronic controls.

## Technical Data:

Speed range:	200 to 3200rpm
Operating modes:	Touch or continuous
Operating temp. range:	+4°C +65°C
Dimensions:	5 x 6.3 x 6.7 in. (13 x 16 x 17 cm.)
Weight:	8.4 lbs. (3.8 kg)
Electrical:	115V AC, 60Hz, 1.5A 230V AC, 50Hz, 0.75A 230V AC, 60Hz, 0.75A
Warranty:	2 years



Unique Multi-Head allows for standard vortexing PLUS holds up to 8 microtubes



## Cat. No.

## Description

ES59382	Mortexer Vortex Mixer With Multi-Head™
ES3370	3" flat head
ES3371	Combo head for microplate and microtubes
ES59953	Horizontal head for vortex mixer (old style), for 12 x 1.5mL
ES59954	Horizontal head for vortex mixer (old style), for 4 x 15mL
ES59955	Horizontal head for vortex mixer (old style), for 2 x 50mL