

AzerPette™ Multi-Channel Pipettors

Ideal for all 96 well applications, the new AzerPette™ Multi-Channel Pipettes are light weight, rugged, and stand up to the demands of everyday use in a busy laboratory environment. Five overlapping volume options are optimized for maximum precision throughout the 1 to 300µL usable range.

Full volume range (1 to 300µL)

- Ten models with overlapping variable volume ranges are available for maximum precision at every conceivable volume up to 300 µL. The micrometer is infinitely adjustable to increments as small as 0.02µL.

Outstanding user comfort

- Each manifold is rotatable and has an integral “dynamic stepped ejection” feature which expels tips two at a time, in undetectable micro-increments, reducing the force required to eject all 8 or 12 tips. The ejector button has also been carefully designed in an ergonomically contoured shape, minimizing pressure on the thumb.

Compatible with standard tips

- All AzerPette™ Multi-Channel Pipettes are compatible with standard 10µL (including extended length), 200µL and 300µL, as well as all standard filter tips. Precision designed pipette shafts/ cones ensure easy loading and a leak-proof fit.

Superior accuracy and reproducibility

- The finely polished high precision piston and chemically resistant, maintenance-free seals provide unsurpassed accuracy. The air-jacketed handle protects the micrometer mechanism from the inaccuracies that can be caused by thermal expansion. This enhances continued pipetting precision even after prolonged use.

Cat. No.	Description	UOM
ES81943-10	AzerPette™ 8- Channel Pipette, 1 to 10 µL	Each
ES81943-20	AzerPette™ 8-Channel Pipette, 2 to 20 µL	Each
ES81943-50	AzerPette™ 8-Channel Pipette, 5 to 50 µL	Each
ES81943-200	AzerPette™ 8-Channel Pipette, 20 to 200 µL	Each
ES81943-300	AzerPette™ 8-Channel Pipette, 30 to 300 µL	Each
ES81944-10	AzerPette™ 12-Channel Pipette, 1 to 10 µL	Each
ES81944-20	AzerPette™ 12-Channel Pipette, 2 to 20 µL	Each
ES81944-50	AzerPette™ 12-Channel Pipette, 5 to 50 µL	Each
ES81944-200	AzerPette™ 12-Channel Pipette, 20 to 200 µL	Each
ES81944-300	AzerPette™ 12-Channel Pipette, 30 to 300 µL	Each

Pipettor Stand! Fits Multi-Channel Pipettors!

- Universal- fits all popular pipettes
- Non-skid rubber feet
- Multi-channel “window” slots

Cat. No.	Description	UOM
ES81998-3	Pipette Rack, for 3 pipettes (or up to one multi-channel)	Each
ES81998-5	Pipette Rack, for 5 pipettes (or up to two multi-channel)	Each
ES81998-6	Pipette Rack, for 6 pipettes (or up to four multi-channel)	Each

