

Pathology Containers

Snap Top Pathology Containers

- Manufactured from translucent, high-density polyethylene and natural polypropylene
- Sturdy, light-weight, and shatter resistant
- Excellent resistance to histology reagents
- Snap-on lid creates a secure seal



Cat. No.	Style	Size	Diam.	Height	Material	UOM
ES3700	Snap Top	4 oz.	3"	2.25"	HDPE	250/case
ES3701	Snap Top	8 oz.	3.5"	2.375"	PP	250/case
ES3702	Snap Top	16 oz.	3.5"	3.5"	PP	100/case
ES3703	Snap Top	32 oz.	4.75"	4.75"	HDPE	100/case
ES3704	Snap Top	64 oz.	6.5"	4.5"	PP	50/case
ES3705	Snap Top	85 oz.	6.5"	6.5"	PP	25/case
ES3706	Snap Top	165 oz.	9"	6.5"	PP	25/case

Screw Top Pathology Containers

- Manufactured from high-quality, clarified polypropylene
- Screw top lid provides a leak-resistant closure to help protect specimen and increase safety of lab personnel
- Excellent resistance to histology reagents



Cat. No.	Style	Size	Diam.	Height	UOM
ES60120	Screw Top	4 oz.	2.73"	2.3"	50/case
ES60180	Screw Top	6 oz.	2.73"	2.74"	50/case
ES60240	Screw Top	8 oz.	2.73"	3.49"	24/case
ES60360	Screw Top	12 oz.	3.54"	3.20"	24/case
ES60480	Screw Top	16 oz.	3.54"	3.85"	24/case
ES61000	Screw Top	32 oz.	4.47"	3.84"	12/case

Large Screw Top Pathology Containers

- Manufactured from high-quality, high-density polyethylene
- Includes heavy-duty screw cap lid
- Moisture and tamper resistant
- Each container has a securely attached handle for carrying

Cat. No.	Style	Size	Diam.	Height	UOM
MM100	Large Screw Top	0.6 gallon	6.5"	7"	12/case
MM200	Large Screw Top	1.25 gallon	8"	8.5"	12/case
MM350	Large Screw Top	3.5 gallon	12"	11.5"	4/case
MM500	Large Screw Top	5 gallon	12"	16.5"	4/case



Small Screw Top Pathology Containers

Our small screw top pathology containers are available in sizes ranging from 20 mL to 120 mL.

- Polypropylene containers with polyethylene screw caps
- Containers are 95 kPa compliant and safe for pneumatic tube systems
- Our 20 mL and 40 mL containers feature our Formalock™ closure; two clicks ensure the seal is fully engaged for safe collection and transport of specimen

Cat. No.	Style	Size	Diam.	Height	UOM
PFNBF-20EML	Screw Top	20 mL	30mm	40mm	200/case
PFNBF-40EML	Screw Top	40 mL	45mm	40mm	200/case
PFNBF-60EML	Screw Top	60 mL	45mm	58mm	200/case
PFNBF-90EML	Screw Top	90 mL	45mm	73mm	200/case
PFNBF-120EML	Screw Top	120 mL	50mm	76mm	200/case

