



FEEL THE SW DIFFERENCE®


CHEMICAL PROTECTION FOR TECHNICAL ENVIRONMENTS.



NEOPRENE POWDER-FREE SINGLE-USE GLOVES

HIGH-PERFORMANCE NEOPRENE GLOVES PROVIDE SUPERIOR RESISTANCE TO CHEMICALS, PUNCTURES AND ABRASION

- Silicone-free for critical environments
- Exceeds 1.5 AQL laboratory standards
- SHA Dermatologically Approved



Latex-free pure neoprene

Optimal thickness for durability and tactile sensitivity

Tear-resistant beaded cuff for easy donning

PC

WEIGHT	LENGTH	THICKNESS	COLOR
4.4g (large size)	9.5"	2.8 mil	GREEN



Feel the difference.™

SW®
33278 Central Avenue, Unit 102, Union City, CA 94587, USA
Tel: +1.510.429.8692 | Fax: +1.510.487.5347

Trademarks and registered trademarks are the property of SW and its affiliates.
©2023 SW. All rights reserved.



WOMEN OWNED™



TECHNICAL DATA



DISTRIBUTED BY:


Azer Scientific

Serving Scientists Who Change The World®

701 Hemlock Road,
Morgantown, PA 19543
Phone: 1 (877) 770-2937
Email: info@azersci.com
AzerSci.com



NEOPRENE POWDER-FREE SINGLE-USE GLOVES



PRODUCT DESCRIPTION

SPECIFICATIONS

PRODUCT	PowerChem®
GRADE	Powder-Free Exam
MATERIAL	Neoprene
THICKNESS FINGER	3.5mil 0.09mm
PALM	2.8mil 0.07mm
CUFF	2.4mil 0.06mm
COLOR	Green
SURFACE	Textured
LENGTH	9.5" 242mm
WEIGHT	4.4g
FIT TYPE	Ambidextrous
CUFF TYPE	Beaded
WATERTIGHT AQL	1.5
ORIGIN COUNTRY	Indonesia

PHYSICAL PROPERTIES

	PRE-AGING	POST-AGING
TENSILE STRENGTH:	14 MPA	14 MPA
ELONGATION:	500%	400%

* Specifications based on size Large glove

STANDARDS & CERTIFICATIONS

EN ISO 21420:2020
EN 455-2000-2015
ASTM D6977-19
REACH EU 1907/2006
EU 2016/425
FDA 21 CFR 177.2600



ISO 374-1:2016/TYPE C



Please visit swssglobal.com/symbols for information on standards used on this document.

ORDERING

ITEM NO.	SIZE	PALM WIDTH	PACKAGING
ES86240XS	XS 5½ – 6	3.0" 76mm	100 gloves per box 10 boxes per case
ES86240S	S 6½ – 7	3.4" 86mm	100 gloves per box 10 boxes per case
ES86240M	M 7½ – 8	3.7" 95mm	100 gloves per box 10 boxes per case
ES86240L	L 8½ – 9	4.2" 106mm	100 gloves per box 10 boxes per case
ES86240XL	XL 9½ – 10	4.5" 115mm	100 gloves per box 10 boxes per case
ES86240XXL	2XL 10½ – 11	4.7" 120mm	100 gloves per box 10 boxes per case

APPLICATIONS

- For Medical Use
- Chemical Handling
- Light Chemical Mixing

Keep out of sunlight. Store in a cool, dry place. Keep away from sources of ozone or ignition.



R_K00040X0323



SW®
33278 Central Avenue, Unit 102, Union City, CA 94587, USA
Tel: +1.510.429.8692 | Fax: +1.510.487.5347

Trademarks and registered trademarks are the property of SW and its affiliates.
©2023 SW. All rights reserved.

