

IN-LAB PRINTING SOLUTIONS

AZERPATH™

YOUR SOLUTION FOR
• EFFICIENT • RELIABLE • DURABLE PRINTING



Primo

CASSETTE PRINTER

Specifications:

Dimensions (L x W x H):	20 x 11.8 x 29 in. (510 x 300 x 737 mm)
Weight:	67 lbs. (30 kg)
Power requirements:	AC 110V, 50/60 Hz, 350 Watts
Laser model:	Ultraviolet laser
Laser power:	1 Watt
Working temperature:	59°F - 86°F (15°C - 30°C)
Work noise:	≤ 65dB
Data interface:	USB 2.0
Cassette printing speed:	3 seconds per cassette (on-demand printing)
Resolution:	2,500 DPI
Number of hoppers:	6
Cassette capacity per hopper:	80 cassettes per hopper
Total cassette capacity:	480 cassettes
Output tray capacity:	6 cassettes



Cat. No.	Description
ES39915	Primo™ Cassette Printer

- Ultraviolet laser (ribbon and cartridge free printing)
- 6 cassette magazines hold up to 480 pcs total, load 6 different color cassettes for efficient workflow
- Precision construction for jam-free cassette printing
- Integrate to LIS/HIS System
- Compatible with hinged cassettes and able to print cassettes with and without lids

Primo

BP CASSETTE PRINTER

Specifications:

Dimensions (L x W x H):	20 x 14 x 17 in. (500 x 350 x 435 mm)
Weight:	66 lbs. (30 kg)
Power requirements:	AC 110V, 50/60 Hz, 350 Watts
Laser model:	Ultraviolet laser
Laser power:	1 Watt
Working temperature:	59°F - 86°F (15°C - 30°C)
Work noise:	≤ 65dB
Data interface:	USB 2.0
Cassette printing speed:	3 seconds per cassette (on-demand printing)
Resolution:	2,500 DPI
Number of hoppers:	6
Number of output trays:	2
Cassette capacity per hopper:	80 cassettes per hopper
Total cassette capacity:	480 cassettes
Output tray capacity:	150 cassettes
Total output tray capacity:	300 cassettes



Cat. No.	Description
ES39916-BP	Primo™ BP Cassette Printer

- Ultraviolet laser (ribbon and cartridge free printing)
- Software system enables high-resolution and high content 1D and 2D barcodes, alphanumeric, and graphic
- On-demand or batch mode printing, high print speed of up to 20 cassettes per minute
- Install scanner for instant scanning
- Convenient to be used at workstation is a portable output tray holds up to 150 cassettes

Cassette Printer Consumables



ES6014




ES6015



ES6016



ES6047



ES6048


Slimsette - Lids Hinged at Front

Tissue 2,000/case	Cat No.	Biopsy 2,000/case	Cat No.	Mesh 2,000/case	Cat No.
Aqua	ES6015-A	Aqua	ES6016-A	Aqua	ES6014-A
Blue	ES6015-B	Blue	ES6016-B	Blue	ES6014-B
Green	ES6015-G	Green	ES6016-G	Green	ES6014-G
Gray	ES6015-GY	Gray	ES6016-GY	Gray	ES6014-GY
Lilac	ES6015-L	Lilac	ES6016-L	Lilac	ES6014-L
Orange	ES6015-O	Orange	ES6016-O	Orange	ES6014-O
Pink	ES6015-P	Pink	ES6016-P	Pink	ES6014-P
Peach	ES6015-PH	Peach	ES6016-PH	Peach	ES6014-PH
Tan	ES6015-T	Tan	ES6016-T	Tan	ES6014-T
White	ES6015-W	White	ES6016-W	White	ES6014-W
Yellow	ES6015-Y	Yellow	ES6016-Y	Yellow	ES6014-Y

Hinged, Separate Lids

Tissue 2,000/case	Cat No.	Biopsy 2,000/case	Cat No.
Aqua	ES6047-A	Aqua	ES6048-A
Blue	ES6047-B	Blue	ES6048-B
Green	ES6047-G	Green	ES6048-G
Gray	ES6047-GY	Gray	ES6048-GY
Lilac	ES6047-L	Lilac	ES6048-L
Orange	ES6047-O	Orange	ES6048-L
Pink	ES6047-P	Pink	ES6048-P
Peach	ES6047-PH	Peach	ES6048-PH
Tan	ES6047-T	Tan	ES6048-T
White	ES6047-W	White	ES6048-W
Yellow	ES6047-Y	Yellow	ES6048-Y

Primo Benefits



EASY INTEGRATION

Simple step-by-step instruction allows you to integrate this printer's software into any LIS/HIS system. The perfect patented-designed technology without the malware to follow in tow.



SWAP HOPPERS EFFORTLESSLY

With a capacity of 80 cassettes per hopper and high speed printing capabilities, the 6-hopper design makes it effortless and fast to swap out. Allows for a variety of cassettes, including with and without lids as well as multiple color changes.

The innovative AzerPath Primo™ Cassette Printers will help to improve both efficiency and productivity in any histology laboratory while reducing technician downtime.

Learn More


Scan here to see all the product listing of our AzerPath Primo™ Cassette Printers or Presto™ Slide Printers. Find more information or schedule demos!

Presto Benefits



INNOVATIVE PRINT QUALITY

With UV laser technology, the AzerPath Presto™ Printers delivers high-resolution (2,500 DPIs) printing that is permanent and resistant to degradation. Print alphanumeric, graphic, logo, QR codes and barcodes.



ELIMINATE COSTLY CONSUMABLES

Thanks to the UV laser technology, you'll be able to say goodbye to those costly printer consumables. No more hassling with print heads or ribbons and print cartridges.

The cutting-edge AzerPath Presto™ Slide Printers provide an effective solution to time management workflow while assuring easy integration in any laboratory or hospital information system.



Presto

UV MICRO SLIDE PRINTER



Specifications:

Dimensions (L x W x H):	16 x 7 x 15 in.(395 x 188 x 372 mm)
Weight:	40lbs. (18kg)
Power requirements:	AC 110V, 50/60 Hz, 350 Watts
Laser model:	Ultraviolet laser
Laser power:	1 Watt
Work noise:	≤ 65dB
Data interface:	USB 2.0
Resolution:	2,500 DPI
Slide printing speed:	< 3 sec. per slide, 20 slides per min.
Slide capacity:	100 slides
Output tray capacity:	10 slides

Cat. No.

Description

ES39911 Presto™ UV Micro Slide Printer

- The downward slide delivery system allows smooth and safe feed of slides
- Print is crisp and permanent, resistant to all chemical solvents during processing
- Versatile printout - Alphanumeric, multiple languages, symbols, 1D & 2D codes, & graphics
- Compact size, ideal to place on the top of the microtome
- Built-in barcode reader, scans barcodes quickly

Presto

DUO UV SLIDE PRINTER



Specifications:

Dimensions (L x W x H):	17 x 9 x 14 in. (425 x 235 x 355mm)
Weight:	46lbs. (21 kg)
Power requirements:	AC, 110V, 50/60 Hz, 300 Watts
Laser model:	Ultraviolet laser
Laser power:	1 Watt
Work noise:	≤ 65dB
Data interface:	USB 2.0
Resolution:	2,500 DPI
Slide printing speed:	< 3 sec. per slide, 20 slides per min.
Slide capacity:	200 slides
Output tray capacity:	40 slides

Cat. No.

Description

ES39912 Presto™ Duo UV Slide Printer

- Dual hopper functionality, with two 100 piece hoppers
- Option of different colors or charged/non charged slides
- Switch seamlessly from one hopper to the other
- Includes setup and 2 pre-installed templates
- On-demand and batch mode printing
- Upward delivery system

Slide Printer Consumables

DigiVU™ Microscope Slides

- Pre-cleaned, smooth, grounded edges
- White glass
- Optimized for compatibility with Azer printers
- Quality tested for optimal results for digital pathology



Color		Standard	Standard (+)	UOM
Blue		90DVSL-B	90DVSL-B+	72 slides/box, 20 boxes/case
Green		90DVSL-G	90DVSL-G+	72 slides/box, 20 boxes/case
Lilac		90DVSL-L	90DVSL-L+	72 slides/box, 20 boxes/case
Pink		90DVSL-P	90DVSL-P+	72 slides/box, 20 boxes/case
White		90DVSL-W	90DVSL-W+	72 slides/box, 20 boxes/case
Yellow		90DVSL-Y	90DVSL-Y+	72 slides/box, 20 boxes/case

Azer-Ink Chemically Resistant Markers

- Ultra-fine tip leaves a permanent mark on slides and cassettes
- Designed specifically for histology use
- Smear-proof, waterproof, and odorless
- Resistant to chemicals such as xylene and alcohol



Cat. No.	Description	UOM
ES-001	Azer-Ink, chemically resistant markers	12/box

We have your workflow covered.

When it comes to histology, Azer Scientific has everything a lab could need. From prefilled formalin to our Azer antibodies for IHC and everything in between, including:

- A comprehensive selection of buffers, fixatives, specimen collection containers, and more
- Embedding equipment and consumables, like our AzerPath Embedding Center and Azerfin Paraffin
- Industry leading Presto™ and Primo™ UV Printers
- High quality DigiVU™ Microscope Slides for optimal standard imaging and digital pathology
- Histotech-preferred iCON™ H&E Stain for consistent, repeatable staining as well as special stain kits/components and IHC antibodies and reagents

Private Labeling Made Easy

Azer Scientific can customize your prefilled formalin vials or specimen collection kits to ensure your lab's collection is streamlined. From custom labels and packaging to specific kit components, we're here to meet your needs.

Why Choose Azer Scientific?

Azer Scientific strives to continue pushing the boundaries of the ever-evolving world of research by offering solutions fueled by innovation. We specialize in producing products designed to make every experiment / study / diagnosis faster and more efficient. As an ISO-certified manufacturer with an ever-expanding inventory, our aim is to cater to the unique needs of every customer by providing them with superior products at an exceptional value.



POWERED BY

