

AzerPro Repeating Pipettor

- Offers 9 tip sizes, 10 volume adjustments and 120 program steps
- **Maintenance-free:** does not require additional resources such as a power supply or batteries
- Easy-moving dosing button and light-weight design allows for convenient dispensing for long dosing series
- Volume adjustments and dosing volumes are shown in a vision panel located at the top of the pipettor
- Specially designed to precisely align with AzerPro Syringe Tips
- Resistant to common laboratory chemicals
- Every device is tested for functionality and precision

NO ▶	1	2	3	4	5	6	7	8	9	10
0.05 ml	1	1.5	2	2.5	3	3.5	4	4.5	5	5.5
0.10 ml	2	3	4	5	6	7	8	9	10	11
0.20 ml	4	6	8	10	12	14	16	18	20	22
0.50 ml	10	15	20	25	30	35	40	45	50	55
1.00 ml	20	30	40	50	60	70	80	90	100	110
1.25 ml	25	37.5	50	62.5	75	87.5	100	112.5	125	137.5
2.50 ml	50	75	100	125	150	175	200	225	250	275
5.00 ml	100	150	200	250	300	350	400	450	500	550
10.00 ml	200	300	400	500	600	700	800	900	1000	1100
12.50 ml	250	375	500	625	750	875	1000	1125	1250	1375
25.00 ml	500	750	1000	1250	1500	1750	2000	2250	2500	2750
50.00 ml	1000	1500	2000	2500	3000	3500	4000	4500	5000	5500
STEPS ▶	48	31	22	18	14	12	10	9	8	7

Tip size

Cat. No.	Description	UOM
ES42020 -0000	AzerPro Repeating Pipettor	Each



Sapphire Syringe Tips

- Designed for the AzerPro Repeating Pipettor
- Available in nine different tip sizes, for dispensing volumes from 1µl to 50 mL
- Provides reliable, and repeatable results for a wide range of liquids including: viscous, infections, and volatile
- Available in **sterile** and **non-sterile options**
- All sterile tips are individually wrapped and certified to be pyrogen, ATP, RNase, and DNase free

Classic Syringe Tips



Size	NON-STERILE		STERILE	
	Cat. No.	Cat. No.	Cat. No.	UOM
0.05 mL	ES40010-0000	ES40010-0008		100/pack
0.5 mL	ES40010-0001	ES40010-0009		100/pack
1.25 mL	ES40010-0002	ES40010-0010		100/pack
2.5 mL	ES40010-0003	ES40010-0011		100/pack
5.0 mL	ES40010-0004	ES40010-0012		100/pack
12.5 mL	ES40010-0005	ES40010-0013		100/pack

Classic Style Syringe Tips (standard and sterile), are suitable for: Ripette®, Ripette®pro, Ripette®genX, eLine™-Dispenser, Minilab 100, Minilab 101, Minilab 201, Distriman®, StepMate, EasyStep, EasyStep electronic, RepeatOne®, Handrop.

Evolution Syringe Tips



Size	NON-STERILE		STERILE	
	Cat. No.	Cat. No.	Cat. No.	UOM
0.1 ml	ES40017-0000	ES40017-0011		100/Pack
0.2 ml	ES40017-0001	ES40017-0012		100/Pack
0.5 ml	ES40017-0003	ES40017-0013		100/Pack
1.0 ml	ES40017-0004	ES40017-0014		100/Pack
2.5 ml	ES40017-0005	ES40017-0015		100/Pack
5.0 ml	ES40017-0006	ES40017-0016		100/Pack
10.0 ml	ES40017-0007	ES40017-0017		100/Pack
25.0 ml	ES40017-0009	ES40017-0019		25/Pack
50.0 ml	ES40017-0010	ES40017-0020		25/Pack

Evolution Syringe Tips (standard and sterile), are suitable for: Ripette®, Ripette®pro, STEPMATE, Minilab 100/101, Multipette® 4780, Multipette® plus, Multipette® M4, HandyStep®, HandyStep® electronic, Distriman



Adapters

Cat. No.	Description	UOM
ES43001-0847	Adapter for 25 ml Tip	Each
ES43001-0848	Adapter for 50 ml Tip	Each

All product names are registered by their respective owners:

HandyStep® - Brand GmbH & Co. KG | bioclean®, Ripette® - Ritter GmbH | eLine™ - Sartorius AG | RepeatOne® - Starlab AG | Distriman®, Repetman® - Gilson Inc. | Multipette® - Eppendorf AG