

Rotating Pipettor Holder

Holds up to Six Standard Single-Channel Pipettes



This pipette carousel is designed to accommodate a set of six single-channel pipettes and keep them within your reach, ready to use. The freely rotating circular rack section has been painstakingly engineered to fit and secure not only the most popular pipette types, but over twenty different pipette brands and the vast majority of all pipette brands on the market today.

With a convenient handle knob at the top, the entire carousel can be easily and safely carried or repositioned, even with a full load of six pipettes. The wide heavyweight base contains a hidden steel stability plate, providing a solid foundation for the carousel and eliminates any possibility of tipping. Compact (only 6.25 inches in diameter), stable and versatile, this carousel is a welcome addition to any research laboratory.

Features:

- Accepts most popular pipettes
- Stable, sturdy base
- Six pipette capacity
- Rotates 360°

Specifications:

Dimensions (W x H): 6.25 x 12in (16 x 30cm)



Cat. No.	Description	UOM
ES3924	Rotating Pipettor Holder for up to Six Standard Single-Channel Pipettes	Each