



SECURESAFE™

SAFETY BLOOD COLLECTION SET



- Winged blood collection needle bonded to a flexible tubing with Luer connector
- SecureSafe™ safety feature allows for one-handed activation to cover the needle immediately after use to help protect against accidental needlestick injuries
- Color-coded wings for easy identification
- 3/4" needle length; 12" tubing



Cat. No.	Description	UOM
ES70115-21G	21G Collection Set, Green	50/pack
ES70115-23G	23G Collection Set, Blue	50/pack
ES70115-25G	25G Collection Set, Yellow	50/pack