

Serological Pipettes

Sapphire 



Our serological pipettes are manufactured using high-grade polystyrene for excellent for clear observation and low retention of liquid on the pipet surface to assure accurate delivery.

Applications include cell and tissue culture, bacteriological, and clinical research. Sterilized by E-beam radiation, our pipettes are DNase/Rnase/Pyrogen free, and non-cytotoxic.

The graduations are calibrated for accurate dispensing to within $\pm 2\%$.

Cat. No.	Description	UOM
ES8030039	5 mL, serological pipette, sterile, individually wrapped	50/bag, 500/pack
ES8030049	10 mL, serological pipette, sterile, individually wrapped	50/bag, 400/pack
ES8030059	25 mL, serological pipette, sterile, individually wrapped	25/bag, 200/pack
ES8030069	50 mL, serological pipette, sterile, individually wrapped	25/bag, 100/pack