

nUVaClean™ | UV Pipette Carousel with Germicidal UV Lamp

The nUVaClean™ is the first and only pipette rack that not only organizes and protects six pipettes, but also destroys up to 99% of unwanted contamination. The nUVaClean™ uses a germicidal UV-C mercury vapor lamp. A high-efficiency reflector/concentrator completely bathes pipette shafts in 254 nm UV light. This exposure destroys unwanted microbiological contaminants and prevents cross-contamination in PCR and other sensitive procedures.

- Holds up to 6 pipettes
- Eliminates cross contamination
- 360° UV exposure
- Compact, 6.5" footprint
- Cleans 6 pipettes in 28 minutes



Cat. No.	Description	UOM
ES81994	nUVaClean™ UV Pipette Carousel with germicidal UV lamp	Each