

xPIPETTE single channel pipette



Patented ejection technology reduces fatigue and makes the xPipette a top choice to improve your lab's ergonomics. Designed with durable and lightweight components to produce a precise and accurate liquid transfer while maintaining hand comfort.

Features:

- Calibrated to ISO 8655 standards with Biotix xTIP4 tips, also fully compatible with Rainin™ LTS™
- Light spring forces produce a smooth and controlled aspiration with perfect balance
- Padded ejector button reduces the dynamic forces and vibrations your hand is exposed to during tip ejection
- Simple volume adjustment by turning the plunger

* Rainin™ and LTS™ are a registered trademarks of Rainin™.

Cat. No.	Description	UOM
ES81949-2	xPIPETTE single channel pipette, 2 µL	Each
ES81949-10	xPIPETTE single channel pipette, 10 µL	Each
ES81949-20	xPIPETTE single channel pipette, 20 µL	Each
ES81949-100	xPIPETTE single channel pipette, 100 µL	Each
ES81949-200	xPIPETTE single channel pipette, 200 µL	Each
ES81949-1M	xPIPETTE single channel pipette, 1000 µL	Each