

DEPC Treated Water



DEPC-treated water is suitable for use with RNA. It is prepared by incubating with 0.1% diethylpyrocarbonate (DEPC) and is then autoclaved to remove the DEPC and sterile filtered.

Features:

- Contains DEPC 0.1 mL/L
- Molecular Grade Water, sterile
- DNase/RNase and pyrogen free

| Cat No. | Size | UOM |
|-----------|---------|------|
| ES58180-1 | 250 mL | Each |
| ES58180-2 | 500 mL | Each |
| ES58180-3 | 1 Liter | Each |