

Picus® 2 Electronic Pipette

The Next Generation of Pipetting

Picus® 2 Electronic Pipette


- Reliable performance and reproducible results
- Easily accessible pipetting modes and customizable programs speed up repetitive work
- Lightweight and ergonomic design reduces muscle strain and provides comfortable pipetting
- Sartorius pipetting mobile app provides configurable protocols and pipette management
- Open connectivity feature provides integration of Picus® 2 pipette into existing lab operating systems
- Service and calibration reminders
- Password protection and distinct administration user profiles protect against unauthorized modifications
- Calibration adjustment to provide accurate results even with the most challenging liquid
- Easy to clean and decontaminate; lower parts of the pipette can be decontaminated by autoclaving (excluding 1,200 µL multichannel models), with UV, or with disinfectants
- Battery typically lasts a full working day and can be fully recharged in an hour

Sartorius Mobile App:

- Sample preparation workflows and pipette management
- Step-by-step workflow guidance and automatic selection of correct settings for Picus® 2
- Pipette updates
- Ever-growing database of configurable protocols



Learn More



Ordering Information: Picus® 2

Order Codes	Channels	Volume Range (µl)	Safe-Cone Standard	Filters Plus	Optifit Tips (µl)	SafetySpace Filter Tips (µl)
ES81953-2-110	1 	0.2-10	-	-	10	10
ES81953-2-1120	1 	5-120	721008	721018	200	120
ES81953-2-1300	1 	10-300	721007	721017	350	300
ES81953-2-11M	1 	50-1,000	721006	721016	1,000	1,000
ES81953-2-15M	1 	100-5,000	721005	721015	5,000	5,000
ES81953-2-110M	1 	500-10,000	721005	721015	10,000	-
ES81953-2-810	8 	0.2-10	-	-	10	10
ES81953-2-8120	8 	5-120	721008	721018	200	120
ES81953-2-81200	8 	10-300	721007	721017	350	300
ES81953-2-8300	8 	50-1,200	721006	721016	1,200	1,200
ES81953-2-1210	12 	0.2-10	-	-	10	10
ES81953-2-12120	12 	5-120	721008	721018	200	120
ES81953-2-12300	12 	10-300	721007	721017	350	300
ES81953-2-121200	12 	50-1,200	721006	721016	1,200	1,200



Pipette Stands and Charging Options

Order Code	Name
ES81953-CS1	Charging stand for 1 Picus® pipette
ES81953-CS4	Charging carousel for 4 Picus® pipettes