

Prefilled Formalin Containers

Azer Scientific's prefilled formalin containers are half-filled with a traditional 10% neutral buffered formalin fixative, used in histological applications for fixation and storage of tissue specimens. Our line of prefilled containers are specifically designed for safe collection and transport of various sized tissue specimens.

- Polypropylene containers with polyethylene screw caps
- Smaller containers (20 mL-120 mL) are 95 kPa compliant
- Inner packaging allows for easy breakdown of the cases for distribution to other labs and doctors offices
- Inner packaging made of chipboard material, safe for operating room use

Formalock™

Smaller sized prefilled formalin containers feature Formalock™ closures.

- Two audible “clicks” to assure the user that the seal is fully engaged
- Eliminates annoyance and hazard of leaking prefills
- Ensures safe transport of specimen



Prefilled Formalin Containers

Cat. No.	Description	UOM
PFNBF-07	7 mL (5 mL fill)	100/tray, 5 tray/case (500)
PFNBF-20	Formalock™, 20 mL (1/2 filled)	32/pack, 4 packs/case (128)
B-PFNBF-20	Formalock™, 20 mL (1/2 filled), bulk	32/pack, 24 packs/case (768)
PFNBF-40	Formalock™, 40 mL (1/2 filled)	24/pack, 2 packs/case (48)
B-PFNBF-40	Formalock™, 40 mL (1/2 filled), bulk	24/pack, 20 packs/case (480)
B-PFNBF-60	Formalock™, 60 mL (1/2 filled), bulk	27/pack, 12 packs/case (324)
PFNBF-60	Formalock™, 60 mL (1/2 filled)	27/pack, 2 packs/case (54)
PFNBF-90	Formalock™, 90 mL (1/2 filled)	12/pack, 4 packs/case (48)
B-PFNBF-90	Formalock™, 90 mL (1/2 filled), bulk	12/pack, 20 packs/case (240)
PFNBF-120	Formalock™, 120 mL (1/2 filled)	24/pack, 2 packs/case (48)
B-PFNBF120	Formalock™, 120 mL (1/2 filled), bulk	24/pack, 8 packs/case (192)

Prefilled Formalin Containers

Our large prefilled formalin containers are designed for collection of large tissue specimens. Each container is half-filled with 10% neutral buffered formalin.

Containers sized 180 mL-1000 mL containers feature a heat seal to prevent leakage during shipment.

We also offer a 0.6 gallon and 1.25 gallon prefilled container. These containers are made from high-density polyethylene and feature a handle for easy transport. The secure screw top closure creates a moisture and tamper resistant seal.

The 3840 mL container is a metal paint can, with attached metal handle.



Cat. No.	Description	UOM
PFNBF-180	Prefilled Formalin Container, 180 mL (1/2 filled)	50/case
PFNBF-240	Prefilled Formalin Container, 240 mL (1/2 filled)	24/case
PFNBF-480	Prefilled Formalin Container, 480 mL (1/2 filled)	24/case
PFNBF-1000	Prefilled Formalin Container, 1000 mL (1/2 filled)	12/case
PFNBF-0.6G	Prefilled Formalin Container, 0.6 gallon (1/2 filled)	Each
PFNBF-1.25G	Prefilled Formalin Container, 1.25 gallon (1/2 filled)	Each
PFNBF-3840	Prefilled Formalin Container, 3840 mL paint can (1/2 filled)	4/case

Customization Options

Don't see what you are looking for here? We are able to customize your prefilled containers with your choice of container size, and label. Contact us to learn more!