

myBath Water Bath Series


Features:

- Removable, hinged lid*
- Quik-CAL™ user calibration function
- Seamless, stainless steel chamber
- Digital accuracy and convenience

Tube racks included with myBath 4L!


The myBath series is designed to guarantee superior accuracy and performance in a compact footprint. With four sizes to choose from, there is a myBath to fit every lab.

The 4, 8, and 12L models include easy-to-use controls that allow the user to digitally select and monitor settings. Adjustable in 0.1°C increments, the temperature is shown on a large, 3-digit display. The proprietary Quick-CAL™ feature allows for simple, in-lab calibration.* High efficiency insulation maintains temperature stability and uniformity within the bath, while keeping the housing cool to the touch.

All myBath models feature seamless, stainless steel chambers for easy cleaning and long life. The chambers resist rust and have no corners to collect dirt or contamination. The larger 8L and 12L units have drains for hassle-free water changes.

The 4L, 8L, and 12L models include a see-through cover that can be removed, or flipped back on a rear hinge to conserve bench space. The lid's gabled design will divert condensation to the front and back of the bath, away from samples.

Racks are available to accommodate 1.5/2.0mL, 15mL, and 50mL tubes within the 4L, 8L and 12L baths. The myBath 4L is supplied with all three racks.

	2 Liter	4 Liter	8 Liter	12 Liter
Cat No.:	ES89551	ES89553	ES89555	ES89557
Temp. Range:	Ambient +5 to 100°C	Ambient +5 to 100°C	Ambient +5 to 100°C	Ambient +5 to 100°C
Temp. Accuracy:	± 0.2°C at 37°C	± 0.2°C at 37°C	± 0.5°C at 37°C	± 0.5°C at 37°C
Temp. Increments:	1.0°C	0.1°C	0.1°C	0.1°C
Temp. Uniformity:	± 0.2°C at 37°C	± 0.2°C at 37°C	± 0.2°C at 37°C	± 0.2°C at 37°C
Exterior (W x D x H):	6.5 x 8 x 6 in. (16.5 x 20.3 x 15.2 cm.)	7.7 x 13.9 x 9 in. (19.5 x 35.3 x 22.9 cm.)	17.5 x 10.5 x 12.4 in. (44.5 x 26.7 x 31.5 cm.)	18 x 10.5 x 16.2 in. (45.7 x 26.7 x 41.1 cm.)
Chamber (W x D x H):	5.5 x 6 x 4 in. (13.9 x 15.2 x 10.2 cm.)	5.3 x 9.4 x 5.5 in. (13.9 x 23.9 x 13.9 cm.)	12.8 x 10.4 x 5.9 in. (32.5 x 26.4 x 14.9 cm.)	32.5 x 10.4 x 7.9 in. (32.5 x 26.4 x 20.1 cm.)
Weight:	4.5 lbs (2.4 kg)	15 lbs (6.8 kg)	17 lbs (7.7 kg)	17 lbs (7.7 kg)
Electrical:	115 or 230V, 50-60Hz, 250W-700W			

Cat. No.	Tube Rack Desc.	4L Rack Capacity	8/12L Rack Capacity
ES89559-R4005	40 x 0.5mL tubes	1	2
ES89559-R40015	48 x 1.5/2.0mL tubes	1 (Included)	2
ES89559-R4115	41 x 15mL tubes	1 (Included)	0
ES89559-R1550	15 x 50mL tubes	1 (Included)	0
ES89559-R7615	76 x 15mL tubes	0	1
ES89559-R3050	30 x 50mL tubes	0	1



myBath 2L

*2L models use traditional calibration and include flat, removable lid.