

Weigh Boats



Made of antistatic polystyrene and will withstand temperatures to 90°C, as well as exposure to diluted acids, alcohols, bases and aqueous solutions. Their flexibility and smooth interior surface with rounded corners allow for easy decanting of powdered samples. Their flat bottoms ensure stability and prevent tipping.

Originally designed as weighing containers, they have proven ideal for dispensing & storage and as disposable mixing/quick freezing trays. They are available in black or white for improved sample contrast/visibility and packaged in ultra-compact dispenser packs.

Features:

- Antistatic and biologically inert
- Smooth surface/rounded corners
- Recessed corners; precise pouring
- Wide, flat bottoms - resist tipping
- For liquids, powders or granular samples
- White or black for visual contrast

Cat. No.	Description	Color	UOM
ES58542-BK	Weigh Boat, 20mL, 45 x 45mm	Black	500/pack
ES58542-W	Weigh Boat, 20mL, 45 x 45mm	White	500/pack
ES58544-BK	Weigh Boat, 100mL, 80 x 80mm	Black	500/pack
ES58544-W	Weigh Boat, 100mL, 80 x 80mm	White	500/pack
ES58546-BK	Weigh Boat, 250mL, 140 x 140mm	Black	500/pack
ES58546-W	Weigh Boat, 250mL, 140 x 140mm	White	500/pack